




GE Medical Systems


gemedicalsystems.com

LOGIQ 200 PRO/α200 SERIES TRANSDUCER GUIDE


H46022LB Linear LB




H46022LH Linear LH




H46022LI Linear LI




H45202LM Linear 10L*




H46022LD Linear LD




H46022LT Linear LT



H46022LE Linear LE



H45202D Sector S317*



Applications	Abdomen OB Gyn	Biopsy Abdomen	Small Parts Vascular	Intraoperative	Intraoperative	Abdomen OB Gyn, Urology	Small Parts Mammography	Cardiac Abdomen
B-Mode Transmit Frequency (MHz)	3.0, 3.5	3.0, 3.5	6.0, 7.5	6.0, 7.0	6.0, 7.0	6.0, 6.5	6.0, 7.0, 9.0	3.0, 3.5, 4.0
FOV	94 mm	94 mm	63 mm*	35 mm	35 mm	83 mm	40 mm	90 degrees
Foot Print (mm)	109 x 22	112 x 267	74 x 16	51 x 20	66 x 21	95 x 21	49 x 15	28 x 22
Biopsy Guide Available		Multi Angle, Reusable	Multi or Single Angle, Reusable				Multi Angle, Reusable	

